

A guide for:

CoJack-Plus Multi-Purpose Table


SHOR-LINE®

▲WARNING This product is for animal use only. We do not condone this product for other uses.

▲WARNING To reduce the risk of injury, the user must read and understand this guide before installation of this product.

{Save This Guide For Future Reference}

PN: 006.3130.01
Updated: 2011



Introduction

Thank you for purchasing Shor-Line products.

We hope this guide answers any questions you might have in regards to use, care or installation. If you can't find the answers to your questions, have further questions or would like to purchase additional products, please contact us at 800.444.1579 or 913.281.1500. If in Europe phone +44 1446 77 20 41.

Consider keeping this guide in a safe and convenient place for future assistance.


IMPORTANT: Mount table to wall completely assembled. DO NOT ATTEMPT TO DISASSEMBLE TABLE WITHOUT CONTACTING SHOR-LINE.

For your safety, please follow the instructions on these labels:

▲CAUTION ▲WARNING

▲WARNING This product is for animal use only. We do not condone this product for other uses.

▲WARNING To reduce the risk of injury, the user must read and understand this guide before installation of this product.

**WARNING** **California Proposition 65 Warning**

This product can expose you to chemicals including nickel (metallic), which are known to the state of California to cause cancer.

For more information go to: www.P65Warnings.ca.gov

{Shor-Line}
511 Osage Ave.
Kansas City, Kansas 66105, USA
phone 800.444.1579 fax 913.281.5339
email guides@shor-line.com
www.shor-line.com

{Shor-Line Limited}
Vale Business Park
Llandow,
Vale of Glamorgan CF71 7PF
United Kingdom
phone +44 1446 77 20 41 fax +44 1446 77 36 68
email quality@shor-line.co.uk
www.shor-line.co.uk

▲CAUTION


Your CoJack-Plus Multi-Purpose Table is heavy so precautions must be taken. Use two or more people when unpacking your shipment.

▲WARNING

Shor-Line products are designed for animal use only. Never allow a person on your CoJack-Plus Multi-Purpose Table.

▲CAUTION

Always set the brakes on your CoJack-Plus Multi-Purpose Table when placing an animal onto the tabletop or when raising or lowering.



{Contents}

Section 1: Unpacking & Inspecting Shipment

Unpacking & Inspecting Steps	4
Table Features	4

Section 2: Operation And Functionality

Adjusting Table Height and Tilt	5
Battery Features	6
Battery Pack Assembly	6
Recharging	6
Control Box	7

Section 3: Battery Pack Specifications

Battery Pack Parts & Assembly	7
-------------------------------------	---

Section 4: Stainless Steel Care & Cleaning

Quality & Life of Your Product	8
Care & Cleaning Best Practices	8
Protecting Against Disease & Infection	8

Section 5: Services, Terms & Conditions, Warranty

Terms & Conditions	9
--------------------------	---



Section 1: Unpacking & Inspecting Shipment

{Unpacking & Inspecting Steps}

Step 1

Locate the Packing Slip that describes the parts and quantities shipped. Your CoJack Plus Multi-Purpose Table should be carefully examined upon arrival for possible damage in shipment. If any damage is found, notify Shor-Line immediately.

Step 2

Review Section 5: Shor-Line's Warranty Information. For more information visit: www.shor-line.com

Step 3

Remove the Table from the shipping carton by cutting the carton away from the table with a utility knife. Due to the weight of the table, we suggest removing the packing material and rolling the table from the carton onto the floor.

Be sure that all parts are contained for your CoJack-Plus Multi-Purpose Table.

Step 4

Check content for concealed damage.


Call Shor-Line immediately so we can expedite replacements or repairs.




Ready to Use


Your CoJack-Plus Multi-Purpose Table arrives ready for use. The battery is fully charged so when you're ready to use your new table simply step on the left foot switch to raise the table and right foot switch to lower.

IMPORTANT: The table will not raise, lower or tilt when the power cord is plugged in for charging.






{Table Features}




The table top will rise to 40" and lower to 14" from floor in 21 seconds with full load (250 lb maximum) enabling you to lift larger breeds with less exposure to back strain and injury.



60° Table Tilt Capability



Step Down
to Set Brake



Section 2: Operation & Functionality

{Adjusting Top Height and Tilt}

Your CoJack-Plus can be used for multi-purposes by adjusting its table top height and tilt.

Step 1

The position of the Rocker Switch determines the functionality of the Foot Switch.

Toggle the switch: (Fig. 2A)

- Left to enable Table Height adjustments
- Right to enable Table Tilt adjustments

Step 2

Adjust Table Height using the Foot Switch.

Step on the:

- Left foot switch to raise the table
- Right foot switch to lower the table

Step 3

Adjust Table Tilt using the Foot Switch. (Fig. 2B)

Step on the:

- Left foot switch to increase the tilt
- Right foot switch to decrease the tilt and return to horizontal position

{Battery Features}

- High battery capacity: 24 Volt DC 2.9 Ah (Approximately 60 lifts at full load 200 lbs.)
- Fast Charge (Approximately 4 hours full charge)
- Medically approved/certified/listed according to EN 6060-1, CANICSA-C22.2 no. 601.1-M90 and UL 2601-1

{Battery Pack Assembly}

See Fig 2C and descriptions below:

1. Control Box Charger
2. Battery Pack
3. Mounting Bracket
4. Battery Latch

NOTE: Item 1 and Item 3 are attached to the table. Item 4 will unlatch to remove Item 2 (Battery Box) from the table.




{Recharge Battery}


Step 1

Activate the foot switch to raise the table so the battery charger is in view.

Step 2

Plug the power cord into the table and the other end into a 120/220 Volt AC outlet. The battery will charge as long as the charge light is on. When the charge light goes out, the battery is fully charged. (Item #6, Fig. 2E)







Step 3


Remove the power cord and the table is ready for use. To maintain a long battery life, plug in charger when table is not in use. (Fig. 2D)

{Control Box Charger} (Fig. 2E)

1. Control Box
2. 120 Volt AC Input
3. Actuator Output 24 Volt DC
4. Foot Switch Control
5. On-Light (only when plugged into 120 Volt AC)
6. Charge Light (only when battery is charging)



Section 3: Battery Pack Specifications



Battery Pack Parts & Assembly			
Item #	Qty.	Part #	Description
1	1	019.7010.00	CONTROL / CHARGER
2	1	070.0900.00	BRACKET MOUNTING MBJ2
3	1	062.0200.05	M5X8 SOC HD CAP SCREW
4	1	159.2400.00	BATTERY PACK WITH BATTERY LATCH



Section 4: Stainless Steel Care & Cleaning

{Quality & Life of Your Product}

IMPORTANT: Routine maintenance will ensure quality and life of your long-term investment and protect animals from transmittable disease and infection

IMPORTANT: To safeguard your investment and rights to the Shor-Line's Warranty, it is important that you know factors that degrade stainless steel. Corrosion, pitting, and discoloration can occur if:

- Bleach, deodorizing agents, disinfectants, and sanitizers (chlorides, bromides, iodides and thiocyanates) stand wet on your product for a prolonged period (over two minutes).
- Bleach, deodorizing agents, disinfectants, and sanitizers are used frequently without being neutralized with clean water after each application.
- Acidic urine is not removed rapidly and surfaces are not washed thoroughly.
- Animal hair, holding urine and fluids, is allowed to collect around legs, hinges, casters, and welded wire intersections.
- Solutions like liquid bleach and antibacterials are allowed to evaporate and dry on surfaces.
- Particles from steel wool, metal fiber brushes, or metal scrapers (which are never recommended for cleaning stainless) get embedded in your product's surface, rust, and give the false impression that the stainless steel is defective.
- Your water source (local or well) is heavy in iron content and water is allowed to dry on your product.
- High levels of pollution (automotive exhaust, aggressive gases, soot, iron oxides, chemicals, etc.) or chloride (from sea or de-icing salts) exist in the environment and frequent, thorough cleaning (manual or rain-washing) does not occur.

{Care & Cleaning Best Practices}


IMPORTANT: You play an essential role in safeguarding your investment and rights to Shor-Line's Warranty. We strongly recommend:

- Daily or weekly cleaning of stainless steel to prolong its longevity and value and to help you maintain its finish.
- Using the mildest cleaning procedure that will do the job effectively.
 - Ordinary waste deposits and fluids can usually be removed with soap and water. Rinse thoroughly with clear water and dry completely with a soft cloth to discourage water spotting.
 - Minor scale build up and some hard water spotting can be removed by washing with vinegar followed by a clearwater rinse and thorough drying.
 - If rubbing, brushing, or scraping is required, use only stainless steel wool, brushes with plastic, nylon or stainless steel fibers, and scrubbers made of plastic or nylon. Ordinary steel wool or steel brushes should never be used.
 - Stubborn deposits or tightly adhering debris may require a commercial cleaning product. Scouring powders will scratch stainless and should not be used. Always review label statements on cleaning products regarding use on metals.
 - Cleaning is not completed until fresh water has been used to thoroughly rinse all residue of salt, bleach, deodorizing agents, disinfectants, and sanitizers.
- Rubbing, drying, and polishing should always be in the direction of the product's grain so scratches will blend into the surface. Efforts opposite the grain may produce a scratched surface appearance.
- Rubbing a soft cloth moistened with olive oil, undiluted white or cider vinegar, or club soda may brighten stainless steel and help remove surface streaks and heat stains.
- Never use paint, lacquer, or varnish as a substitute or shortcut for properly caring and cleaning stainless steel.
- Never intentionally stamp or etch stainless surfaces.

{Protecting Against Disease & Infection}

IMPORTANT: Shor-Line shares your commitment to health and safety. This information is intended help protect your animals, as well as, safeguard your investment.

- Bleaches can be harmful to animals. Bleaches can react with acids and release chlorine gas that is highly toxic.
- Because of its strong corrosive and toxic properties, bleach is diluted with water. The solution's potency (shelf-life) quickly degrades as the active chlorine dissipates with water. Animal care professionals, therefore, cannot know with certainty how effective the bleach solution is over time without implementing a formal monitoring system.
- If chemical deodorizers, disinfectants, or sanitizers are used, protect your animals from respiratory problems and



infections by using clean water to wash residue from all surfaces and dry thoroughly.

- Clorox® Company warns: “We cannot recommend Clorox® to clean metal surfaces, including stainless steel. Prolonged or repeated exposure of metal to strong solutions of our product eventually can cause discoloration or corrosion. To reduce the possibility of metal discoloration, we suggest that exposure time to Clorox® be less than five minutes and that surfaces be rinsed well after application of the solution.”

Section 5: Services, Terms & Conditions

<https://www.midmark.com/service-support/terms-conditions>



SHOR-LINE®

<i>toll free</i>	<i>local</i>	<i>fax</i>
800	913	913
444.1579	281.1500	281.5339

Shor-Line
511 Osage Avenue
Kansas City, KS 66105 U.S.A.

www.shor-line.com