

ANESTHESIA MACHINES

Shor-Line's New Line of Anesthesia

Basic Mobile Anesthesia Machine

Features:

- Five leg base with shielded casters.
- CO2 absorber with pressure gauge.
- "One seal" canister with lifetime warranty against breakage.
- Dual O2 supply ready with check values.
- 0-4 lpm oxygen flowmeter with mechanical stop.
- Oxygen flush with 30 lpm safety restrictor.
- Common gas outlet for non-rebreathing systems.
- Safety pressure relief valve.
- Syringe holder for inflating E/T cuffs.
- Adult & pediatric unilimb rebreathing circuits.
- Non-rebreathing system with 1 liter bag.
- 0.5, 1, 2 and 3 liter bags.
- Ten-year warranty on parts.
- Vaporizer, Dual "E" Tank Manifold and bottles sold separately.



195.1178.00 Basic Mobile Anesthesia Machine 70 lbs.

Premier and Plus Mobile Anesthesia Units

Features:

- "One seal" canister with lifetime warranty against breakage.
- Sturdy A-frame cart with 18" square footprint.
- Plus model has durable, one-level shelf. Premier model has larger shelf with back bar that contains an organizer for supplies
- CO2 absorber with pressure gauge.
- Dual O2 supply ready with check values.
- 0-4 lpm oxygen flowmeter with mechanical stop.
- Oxygen flush with 30 lpm safety restrictor.
- Common gas outlet for non-rebreathing systems.
- Safety pressure relief valve.
- Syringe holder for inflating E/T cuffs.
- Adult & pediatric unilimb rebreathing circuits.
- Non-rebreathing system with 1 liter bag.
- 0.5, 1, 2 and 3 liter bags.
- Ten-year warranty on parts.
- Vaporizer, Dual "E" Tank Manifold and bottles sold separately.



Premier

195.1179.00
195.1180.00



Plus

Plus 70 lbs.
Premier 70 lbs.



Organizer included on
back of Premier unit.