

FELINE SCALE

Quickly receive an accurate weight as felines snuggle in

Benefits:

This compact scale easily fits on almost any counter and delivers an accurate weight on felines and small pets. Its high sides allow the animal to crouch down, making weighing easier.

Features:

- Push button zero/tare weight.
- Pound/kilogram switch.
- 44 lbs./20 kg. capacity, weighs in 2.5 lbs./2 kg increments.
- Tray removes easily for cleaning.
- Operates with 9V battery, making it easily portable.
- Easy to set up with intuitive functioning.
- **SPECS:** 9"W x 21"L x 2"H lip.
- 905.0400.05



Works for small dogs too!

