

MOBILE ANIMAL LIFT TABLE (MALT)

Mobile Animal Lift Table

Benefits:

Shor-Line's multi-functional Mobile Animal Lift Table gives any exam or treatment that much needed lift. This table also comes available with the K9-W8 scale attached, letting you conveniently weigh patients without moving them.

Features:

- Foot switch engages a quiet, battery-operated rechargeable actuator system that can lift up to 200 lbs.
- Table top raises from 10³/₄" to 37³/₈" high.
- Stainless Steel exam top is 44¹/₂"L by 20"W, with raised edges to prevent fluid run-off.
- Features two tie-down restraint cleats on each side to help restrain restless patients.
- Corrosion resistant reinforced nylon locking casters and rotating non-marring bumpers.



903.3120.00 Table
903.3120.02 Table with Scale

187 lbs.
278 lbs.