

CONTINUUM SHORT V-TOP SURGERY TABLE

Shor-Line's Newest V-Top Offering

Continuum Short V-Top Surgery Table

Benefits:

The new 46" L V-Top has all the features you expect from the standard V-Top with a shorter length. You can position the patient from horizontal to near-vertical without external devices. The panel design provides ergonomic ease of use and ready access to trough cleaning. Gives the dependable operation you expect from Shor-Line. Lifts up to 250 lbs.

Features:

- Two polymer five lobe ergonomic knobs allow the panels to open to 270° for easy access to remove surgical waste. Just clean the trough and reinsert the tray for the next surgical procedure.
- Rail attachment on each side.
- A 3" stainless steel trough permits fluid runoff into a drain outlet.
- Continuum Surgery table base is a distinctive stainless steel hemmed column with a clean and sophisticated design.
- The round, tubular "X" base is designed with enough clearance for feet and other surgical tools that need to be in close proximity of the table.
- The base is powder-coated gun metal grey and will stand up to cleaning and disinfecting agents.
- The Hydraulic Jack allows you to place the table anywhere you want with no cords.
- The hydraulic mechanism lifts the table from 30½" to 41".
- **Specs:** 19"W x 46"L
- **Need it taller?** Optional height extenders available.



903.4800.01 Hydraulic Base

227 lbs.