

CLASSIC SHORT V-TOP SURGERY TABLE

Shor-Line's Newest V-Top Offering

Classic Short V-Top Surgery Table

Benefits:

The new 46" L V-Top has all the features you expect from the standard V-Top with a shorter length. The V-Top acts as a helper during surgery, keeping the patient in the correct position. Check out the extra features in this premium table. Gives you flexibility to go from straight horizontal and level to V-Shaped, near vertical positioning without external devices. You choose how you want to operate.

Features:

- High-quality, durable stainless steel table surfaces.
- Two polymer five lobe ergonomic knobs allow the panels to open to 270° for easy access to remove surgical waste. Just clean the trough and reinsert the tray for the next surgical procedure.
- Two-panel system that adjusts from horizontal to near vertical to hold the patient in the desired position during surgery.
- Two sliding, quick-release, cam cleat tie-down attachments provided on each side of the attachment rail.
- 3" stainless steel trough permits fluid runoff into a drain outlet.
- Table can be lifted from 30½" to 41" via hydraulic mechanism.
- Top can rotate 360 degrees or be stationary.
- Solid stainless steel base that lasts for decades.
- **Specs:** 19"W x 46"L
- **Need it taller?** Optional height extenders available.



903.4800.00 Rotational

227 lbs.