

TUB TABLES

Tub Tables

Benefits:

Our multi-purpose tub tables add both efficiency and user comfort to many of your daily procedures. The customizable storage space allows tailoring of your table to your needs. If you're tight on space the 46" table top fits conveniently into most floor plans, but if extra working room is what you're looking for the 60" models deliver. Add the Stainless Steel Exam Top to your purchase and get the perfect surgery and exam table all in one package.

Features:

- Tub tables with knee space come in laminate or Stainless Steel.
- Choose between two storage options: Door and drawers or two doors.
- Standard laminate colors are white or almond. Optional Wilson Art® or Formica® colors are available at additional cost and longer lead time.
- Standard tub tables feature recessed grill covers, and 16" or 6" deep tubs.
- Laminate cabinets are constructed of ¾" MRB (moisture resistant board). door edges are wrapped with 3mm PVC to prevent chipping and add extra protection.
- All tables are 25" wide and 36" high.



904.4000.51



904.5000.09



904.2100.10



904.4000.11

904.2100.10
904.5000.09
904.4000.51
904.4000.11
804.3104.03
804.3004.03

46" Long, 6" Deep Tub, White
46" Long Tub Table, Stainless Steel
46" Long Tub Table, Almond
60" Long Tub Table, White
46" Long Exam Top
60" Long Exam Top

804.3004.03

