


HOW TO BUY ANESTHESIA

WHAT TYPE OF UNIT DO YOU WANT?




- Anesthesia Machine + Oxygen Concentrator
You might consider the Pureline M6000, which is a complete unit with everything but the vaporizer and back-up oxygen.
- Anesthesia Machine
We offer choices for wall mount, table top and mobile.
- Oxygen Concentrator
Are you tired of paying every month for bottles? It may be time to invest in an O₂ concentrator.

MOUNT




- Wall Mount
If you have O₂ connections on the wall, this could be your best choice. We offer both flush mount and mounted with articulating arm.
- Table Top
This gives you flexibility to move to different stations but you will need space to place the unit. (This is most often used in laboratories.)
- Mobile
There are several choices, depending on how you want to provide your oxygen.

OXYGEN



- Oxygen Concentrator
- E-Tank Manifold
(Note: Shor-Line does not sell oxygen)

VAPORIZER



- Isoflurane, Pour or Key Fill
- Sevoflurane, Pour or Key Fill