

FEAR FREE EXAM TABLE

The New Era of Pet Comfort

Features:

- Developed in conjunction with Dr. Marty Becker and Fear FreeSM, which aims to provide a relaxing environment to encourage calmer, more accepting patients.
- Heated Stainless Steel exam top provides a warm surface that has safety shutoffs.
- Non-slip stairs (four steps) fold down gently to the floor with rubber feet, decreasing noise. It operates with fingertip ease, allowing you to focus on the pet.
- 4½" W clips at either end of the table let you attach a washable towel in a soothing color.
- Toe kicks (3"H) around the sides and front of table allow you to comfortably be closer to the table.
- Durable Shor-Line craftsmanship and easy-to-clean materials provide many years of use.
- The Fear FreeSM badge indicates this table is part of a Fear FreeSM practice.
- Comes standard in Light Blue. Other colors available on request with longer lead time.
- **SPECS:** 27"W x 40"L x 36"H.

Benefits:

This new stationary exam table takes the next step in giving the pet a gentle, calming exam experience. The dog can choose to walk up secure, fold down stairs to a warm exam top with a comfy towel. Position the table in peninsula style so you and the pet parent can praise the pet while discussing its health.



903.1800.10 Fear Free, Light Blue

