

MARK III LS EXAM TABLE

This heated top lets you give your patients a warm welcome at exam time.

Mark III LS Exam Table

Benefits:

This exam table is part of our low stress line, combining the versatility of storage with the opportunity to provide a warmer, more inviting exam table experience. The Stainless Steel top has our surgery table's dependable heating technology, keeping the surface an even temperature.

Features:

- Durable Stainless Steel top is formed to encourage drainage to an edge trough.
- Heating unit controller attaches to the side of the table for easy access.
- Controller maintains a factory-set temperature of approximately 93°(33.8°C) with three timed cycle options (60, 90, 120 minutes). (The temperature is NOT adjustable.)
- There is a built-in alarm and shut off to alert you to an overheating situation.
- The cabinet is available in high pressure laminate.
- All cabinets have three drawers and one door.
- Standard laminate color options are White or Almond. Other Wilson Art or Formica colors are available with longer lead times and additional cost.
- Chip-resistant, 3mm PVC edging finishes the laminate door and drawer edges.
- SPECS: Top is 44½"W x 20"D. Storage area is 35.875"W with about a 2" toe kick around the unit. It is 32.875"H.



903.1001.10	Mark III LS Exam Table, White	298 lbs.
905.1001.20	Mark III LS Exam Table, Special Color	298 lbs.
905.1001.50	Mark III LS Exam Table, Almond	298 lbs.

Low-Stress
SHOR-LINE
SINCE 1927