

DRYER CAGES

Shor-Line Grooming Luxury

Dryer Cages

Benefits:

Patented heat-free drying method eliminates burns and dehydration while providing a speedy drying time.

Features:

- Designed with two high volume fans placed on opposite ends of the assembly. These fans act to create turbulence, which lifts and dries the pet's coat while removing the loose undercoat.
- Fan guards at both ends protect paws from fan blades, and multi
- speed settings let you determine how fast the airflow dries patients.
- All fans are UL listed and feature OSHA compliant guards.
- Rounded corners throughout prevent bacterial buildup, making a healthier environment for animals and staff.
- Fluid-retaining lip prevents urine or other leakage.
- Removable PVC coated floors included.
- Wide doors mean easier animal entrance and exit.
- Double unit includes a mobile double large dryer cage on the bottom with another double large dryer cage on the top.
- Includes a removable cage divider gives you the option of drying two smaller animals at one time.
- Triple unit is equipped with stacking plastic cages for holding areas leaving you free to dry the next animal.
- Locking caster base provides instant portability.

SPECS:

- Single Unit is 48" x 36".
- Double Unit is two 48" x 36" cages stacked.
- Triple Unit is one 48" x 36" cage with two smaller holding cages stacked above.
- Each cage compartment is a one-piece mold of high-density polyethylene with a 3/16" thick wall.
- 25½"H x 44¼"W x 28"D. Maximum total air volume is 3,560 CFM.
- Requires 120V AC, 1.4 AMPS per fan.



Tripple Unit



Double Unit



Single Unit

926.4836.01

Single Unit

129 lbs.

926.4836.02

Double Unit

258 lbs.

926.4836.03

Tripple Unit

205 lbs.