

# CAGE DIVIDERS

## Cage Dividers

### Benefits:

Any of our double door cages can instantly be made into two smaller cages with a cage divider. Stainless Steel cages have the option of the Stainless Steel Divider or Thermal Formed ABS Divider.

### Features:

- The Stainless Steel Divider provides heavy-gauge construction.
- Stainless Steel Cages have a removable post system, making it easier to quickly slide recombinant patients into the cage.
- Thermal Formed ABS Divider provides lightweight durability and cleaning ease.
- Thermal Formed ABS handle makes one-handed opening and closing easy.

Removable post



802.0024.00	24"H Cage SS Divider	29"W x 24"H
802.0030.00	30"H Cage SS Divider	29"W x 27"H
802.0036.00	36"H Cage SS Divider	29"W x 33"H
802.0030.01	30"H Cage ABS Divider	30½"W x 27"H
802.0036.01	36"H Cage ABS Divider	30.688 W x 33"H