

MULTI-PURPOSE ULTRASOUND TABLE

Turn all your Shor-Line tables into Ultrasound tables

Benefits:

The Multi-Purpose Ultrasound Table provides flexibility for busy practices to administer ultrasound exams where they are needed. The sturdy, durable table ergonomically supports the pet and the operator, and it folds up easily for storage or to use in another room.



Features:

- Crafted from a lightweight, water-resistant substrate, we include a ½”-thick foam cushion and cover it with a strong polyurethane top. Durability meets comfort for patients and operators.
- The polyurethane top stands up to standard veterinary cleaning procedures.
- Fits standard veterinary table sizes (46”L or 60”L).
- Provides built-in flexibility with the option of moving the legs closer to the edge.
- Combine the small and the large tables to manage large dogs with three opening options. (The legs are positioned 2” inside the top and bottom ends for stability when scanning a large dog on a 60”L table.)
- By adjusting the legs, you can achieve up to a 2” tilt.
- Oval opening allows easy scanning. The large unit has two 6”W x 10½”L openings and the small has one 4”W x 6”L opening.
- Stainless Steel hardware is durable and stands up to cleaning.
- It is available in Small or Large to fit most pets. The large table is intended for medium to large canines and the small is sized for felines and small pets.
- Designed for animal use only.

903.3960.00

903.3961.00

903.3962.00

Small

Large

Both



Specs:

- **Small:** 14”W x 24”L x 11½”H (when legs are fully extended), weighs 9.15 lbs.
- **Large:** 48”W x 24”L x 11½”H (when legs are fully extended), weighs 22.35 lbs.
- **Combined:** 62”W x 24”L x 11½”H (when legs are fully extended) weighs 31.5 lbs.