

FELINE LS SURGERY/EXAM TABLE

Elegance meets versatility in our newest exam table

Benefits:

The low stress Feline LS Surgery/Exam Table creates a warm, welcoming atmosphere for feline exams. The large oval table has plenty of space to set the cat carrier down as well as your laptop or chart.

Features:

- Stylish table has about 120 square inches of space for exams.
- The top has an attractive, rounded design so no one is hitting sharp edges.
- There are two table options:
 - The Blue-Line top has radiused corners and laminated polyurethane edging that holds up to a liter of fluid. It comes in Standard Blue, Jade, Taupe Gray or Black.
 - The T-Mold edge has polyurethane edge molding for a smooth, attractive finish.
- The “X” base provides good room for both you and the pet parent to be close to the table, offering comfort.
- The table jack smoothly adjusts the height from 29” to 41” H so each exam can be at the ideal height.
- The table top is durable, chemical-resistant Stainless Steel that is easy to clean.
- Shor-Line tables are designed for animal use only.

Specs:

- The top is oval-shaped and 30”W x 40”L..
- The hydraulic table base has a jack that lifts the table from 29” to 41”H.
- The based comes with adjustable leveling feet that helps keep the table level on stone or other uneven floors.
- Lifting capacity is 100 lbs.
- Product weight is 225 lbs.

Table Top Options

T-Mold Top

903.4123.00	Standard Blue	T-Mold Edge
903.4123.01	Black	T-Mold Edge

Blue-Line Top

903.4122.00	Standard Blue	Molded Top
903.4122.01	Jade	Molded Top
903.4122.02	Taupe Gray	Molded Top
903.4122.03	Black	Molded Top



903.4123.01

Black

T-Mold Edge



903.4122.00

Standard Blue

Molded Top