

BLUE-LINE SCALES


Shor-Line's Blue-Line Scale

We've made the Blue-Line Scale you love even better while also enhancing your decor. The scales come with new load cell technology and a new combination of materials that delivers incredible accuracy. These products have the ability to withstand years of use.

Blue-Line Scale, Display with Operation Keys

Features:

- Low-profile Stainless Steel platform scale gives accurate and repeatable weights in pounds and/or kilograms.
- Unique polyurethane edge available in four colors.
- Large, easy to see backlit LCD screen with $\frac{3}{4}$ " H numbers.
- Quick, on-screen weight calibration.
- Wall Mount mounts above the scale or another room.
- Specs: 44 $\frac{1}{2}$ "L x 22 $\frac{1}{2}$ "W x 3"H.
- 10' data cord extension available; includes 6' power cord.



Post Mount Display

Blue-Line Scale, Display

Features:

- Remote display works best when there is a clear sightline or path from the platform unit.
- Display with operation keys works like other K9-W8 scales and the remote displays the animal's weight.



Remote Display

Product	Display	Remote Display
Standard Blue, Wall Mount	905.3010.05	905.3010.08
Standard Blue, Post Mount	905.3014.04	905.3014.06
Jade, Wall Mount	905.3011.05	905.3011.08
Jade, Post Mount	905.3015.04	905.3015.06
Taupe Grey, Wall Mount	905.3012.05	905.3012.08
Taupe Grey, Post Mount	905.3016.04	905.3016.06
Black, Wall Mount	905.3013.05	905.3013.08
Black, Post Mount	905.3017.04	905.3017.06