

DRYER CAGES

The Dryer Cages feature durable housing with a heat-free dryer. The special fans lift up the fur to both loosen the wet undercoat and give a fluffy dry for a professional finish.

Benefits:

Patented heat-free drying method provides protection from burns and dehydration while providing a speedy drying time. US Patent #6,058,886.

Features:

- Designed with two high volume fans placed on opposite ends of the assembly. These fans act to create turbulence, which lifts and dries the pet's coat while removing the loose undercoat.
- Fan guards at both ends protect paws from fan blades, and multi-speed settings let you determine how fast the airflow dries patients.
- All fans are UL listed and feature OSHA compliant guards.
- Rounded corners throughout prevent bacterial buildup, making a healthier environment for animals and staff.
- Fluid-retaining lip prevents urine or other leakage.
- Removable PVC coated floors included.
- Wide doors mean easier animal entrance and exit.
- Can also be equipped with stacking plastic cages for holding areas—leaving you free to dry the next animal.
- Includes a removable cage divider gives you the option of drying two smaller animals at one time.
- Double unit includes a mobile double large dryer cage on the bottom with another double large cage dryer on top.
- Locking caster base provides instant portability.
- **Specs:** Single Unit is 48" x 36". Double Unit is two 48" x 36" cages stacked. Triple Unit is one 48" x 36" cage with two smaller holding cages stacked above. Each cage compartment is a one-piece mold of high-density polyethylene with a 3/16" thick wall. 25½"H x 44¼"W x 28"D. Maximum total air volume is 3,560 CFM. Requires 120V AC. 1.4 AMPS per fan.

Heat-Free for a Fluffy Finish



Double Unit



Triple Unit



Single Unit

Single Unit	926.4836.01	129 lbs.
Double Unit	926.2000.02	129 lbs.
Triple Unit	926.3000.03	129 lbs.