

# PVC COATED CAGE FLOORS

## Benefits:

Keep your patients dry and clean with this removable cage floor.

## Features:

- 14-gauge perforated sheet.
- Openings are approximately 3/8" diameter after PVC coating.
- Non-toxic PVC coating is approximately 1/8" thick and seals entire floor.
- Floor height is approximately 2-1/4".



Product	Part No.	Weight
<b>Floor, 18"W Cage</b>	317.0098.00	11 lbs.
<b>Floor, 24"W Cage</b>	317.0098.01	20 lbs.
<b>Floor, 30"W Cage</b>	317.0098.02	25 lbs.
<b>Floor, 36"W Cage</b>	317.0098.03	28 lbs.
<b>Floor, 48"W Single Door Cage</b>	317.0098.04	33 lbs.
<b>Floor, 48"W Double Door Cage*</b>	317.0098.05	41 lbs.
<b>Floor, 60"W Double Door Cage*</b>	317.0098.45	51 lbs.
<b>Floor, 72"W Double Door Cage*</b>	317.0098.46	57 lbs.

\*These cages require 2 pieces. Part numbers above are 1 piece.

