

PRELUDE LED EXAM LIGHTS

Shor-Line's New Line of Exam Lights

Our development team has listened to your needs and developed cutting edge exam lighting to meet them. The Prelude LED Exam Light gives the flexibility to light up issues during an exam while keeping the room cool and saving money.

Prelude Series LED Exam Lights


Benefits:

The 40,000 LUX light feature a color temperature of 4500° Kelvin for visual tissue clarity and at one-meter working distance features a 6-inch field of illumination. The lights use 11 white light LEDs with individual parabolic reflectors with a bulb life of up to 75,000 hours. Keeps the exam room cooler than halogen lights.

An easy-to-use switch allows for low (50%) or high (100%) settings.

Features:

- High efficiency. In general, LEDs consume only 20% of the power halogen bulbs require.
- Very low heat—cool to work under.
- No replacement bulbs or fuses necessary.
- State-of-the-art spring arms facilitate easy, drift-free positioning to give you maximum flexibility in surgery.
- Low maintenance.
- Product has a 3-year warranty.
- Auto voltage sensing of 120V-240V.



913.7000.00 Single Ceiling Mount
913.7000.01 Single Wall Mount

60 lbs.
36 lbs.