

# K9-W8 + BLUE-LINE SCALES


Two great scales for accurate weights

We've made the K9-W8 scale your reliable partner and jazzed up the Blue-Line scale with color. Both scales come with new load cell technology--four load cells accurately weigh four-pawed patients.

## Blue-Line Scale, Display with Operation Keys

### Features:

- Low-profile Stainless Steel platform scale gives accurate and repeatable weights in pounds and/or kilograms.
- Unique polyurethane edge available in four colors.
- Large, easy to see backlit LCD screen with ¾" H numbers.
- Quick, on-screen weight calibration.
- Wall Mount mounts above the scale or another room.
- Specs: 44 ½"L x 22 ½"W x 3"H.
- 10' data cord extension available; includes 6' power cord.



Product	Display
Standard Blue, Wall Mount	905.3010.05
Standard Blue, Post Mount	905.3014.04
Jade, Wall Mount	905.3011.05
Jade, Post Mount	905.3015.04
Taupe Gray, Wall Mount	905.3012.05
Taupe Gray, Post Mount	905.3016.04
Black, Wall Mount	905.3013.05
Black, Post Mount	905.3017.04

## K9-W8 Scale, Display with Operation Keys

### Features:

- Four ultra low profile load cells at each corner provide accurate and repeatable weights in pounds and or kilograms.
- Large, easy to see backlit LCD screen with ¾" high numbers.
- Tare, recall, lbs/kg, light and free/lock features and quick, on-screen weight calibration.
- Gray, ribbed, non-skid mat provides solid patient footing.
- Post display is powder-coated dark gray and stands 37" high. Slide the slim base under the scale or move it within 4" of the scale.
- Specs: 22 ½"W x 44 ½"L x 3"H.
- Comes with 6' power cord. 10' data cord extension available.



905.5010.43	Wall Mount	53 lbs.
905.5010.47	Post Mount	67 lbs.

# FELINE SCALE

*Quickly receive an accurate weight as felines snuggle in*

## **Benefits:**

This compact scale easily fits on almost any counter and delivers an accurate weight on felines and small pets. Its high sides allow the animal to crouch down, making weighing easier.

## **Features:**

- Push button zero/tare weight.
- Pound/kilogram switch.
- 44 lbs./20 kg. capacity, weighs in 2.5 lbs./2 kg increments.
- Tray removes easily for cleaning.
- Operates with 9V battery, making it easily portable.
- Easy to set up with intuitive functioning.
- **SPECS:** 9"W x 21"L x 2"H lip.
- 905.0400.05



Works for small dogs too!

