

K9-W8 SCALES

Four load cells for accurate weight

We've made the K9-W8 Scale you love better with new load cell technology and materials that deliver incredible accuracy and durability.

K9-W8 Scale, Display with Operation Keys

Features:

- Four ultra low profile load cells at each corner provide accurate and repeatable weights in pounds and or kilograms.
- Large, easy to see backlit LCD screen with $\frac{3}{4}$ " high numbers.
- Tare, recall, lbs/kg, light and free/lock features and quick, on-screen weight calibration.
- Gray, ribbed, non-skid mat provides solid patient footing.
- Post display is powder-coated dark gray and stands 37" high. Slide the slim base under the scale or move it within 4" of the scale.
- Specs: 22 $\frac{1}{2}$ "W x 44 $\frac{1}{2}$ "L x 3"H.
- Comes with 6' power cord. 10' data cord extension available.



Post Mount



Wall Mount

905.5010.43
905.5010.47

Wall Mount
Post Mount

53 lbs.
67 lbs.