

OC8000 SERIES OXYGEN CONCENTRATORS

The revolutionary Pureline OC8000 Series takes oxygen concentration to a new level in the veterinary market.

PURELINE® OC8000 SERIES OXYGEN CONCENTRATORS

Benefits:

The new Pureline OC8000 Series Oxygen Concentrators can deliver on demand up to 8 LPM, allowing it to handle the higher flow rates sometimes needed for oxygen therapy. This makes it ideal for use with oxygen therapy cages. It can be an extra oxygen supply for your anesthesia system as well.

Features:

Its state-of-the-art manufacturing in the USA uses high-quality components. It is built with safety in mind with oxygen purity and flow safety alarms. It is UL tested.

- For veterinary use only
- Intended for both oxygen therapy and anesthesia machines
- Extracts medical grade oxygen (USP 93%) from the air and delivers it up to 8 LPM
- Connects with industry standard oxygen DISS fittings
- Includes 10' high pressure green oxygen hose
- Quiet 52-dBA operation
- Energy-efficient; costs average of 20 cents per day



195.1363.00	OC8000	35 lbs.
195.1364.00	OC8200	70 lbs.